

ERCPC

Procedure Table



Marketed & Promoted by

Innovative Therapeutics Pvt Ltd

8/467/2A1, Nehru Street, 2nd Cross Road, 4th Lane, Arthi Nagar,
Saraswathi Nagar, Thirumullaivaoyal, Chennai, Tamil Nadu 600062

+91 44 3559 4099 | +91 93811 27879

sales@itplmedical.com | www.itplmedical.com

Manufactured by

Overview

- HASSLE-FREE ERCP PROCEDURE with our cutting edge ERCP table
- Free C-Arm Movement with height adjustable

Specification

Overall Size	2050 x 850mm (Out Side Size: 2050 X 850 mm)
Patient Surface	1850 x 700mm (800 / 850mm wide option)
Electric Height Range	570 - 940mm.
Hydraulic Height Range	630 - 960mm.
SWL	300kg.
Castor Size	200mm castors and Ezi-Trak 5th wheel steering.
Includes:	Full length side rails that are fully welded stainless steel, emergency CPR backrest release, X-ray sub top option.

From Work

- X-Ray/ C-Arm Comparable 10mm Top
- 32. OD 4mm Pipes /Lifting Supports
40.mmX10.mm Flats
- X-Ray/ C-Arm Comparable 10mm Top

Mattress Work

- Hi Density 3" Imported Foam Mattress,
and Pillow

