

EndoCabi

Endoscope storage &
Drying Cabinet

DQUBE 
An ITPL Group Company



Marketed & Promoted by

Innovative Therapeutics Pvt Ltd

8/467/2A1, Nehru Street, 2nd Cross Road, 4th Lane, Arthi Nagar,
Saraswathi Nagar, Thirumullaivaoyal, Chennai, Tamil Nadu 600062

+91 44 3559 4099 | +91 93811 27879

sales@itplmedical.com | www.itplmedical.com


Manufactured by

DQUBE 
An ITPL Group Company


- **Thoroughly Dry:** Keep your storage environment aseptic and dry through internal air pressure with an average of an air exchange every 20 seconds. Help Prevent microbial growth with channel purging air flow thoroughly drying internal lumens.
- **Vertical is Best:** Prevent scopes from touching and damaging each other, Maintain ideal shape without causing potentially damaging coil memory and Eliminate the chance of water draining onto other scopes.




3 Stage Filtration



20µ Filtration




UV Filtration



10µ Filtration

Technical Specification

Applicable Scopes	Flexible
Capacity	9 Scopes
Air Filtration	3 Stage Filtration
Channel Purging Connectors	Optional
Dimensions	2600 x 620 x 620 mm
Weight	180 Kgs
Inner Wall Construction	Antibacterial HPL Sheet
Outer Wall Construction	Antibacterial HPL Sheet
Doors	Tempered Glass
Lock Type	Manual Lock or Smart Password, Electromagnetic Door Lock
Compressed Air Pressure	2-8 Bar - Optional
Compressed Air Quality	Filtered - Optional
Electrical Requirements	220/230/240 V ~ ; 50/60 Hz



Channel Purging



Designated Slot



HEPA Filter