

# DR-02 EndoWasher

Reliable Automatic Washer disinfector



## Marketed & Promoted by

Innovative Therapeutics Pvt Ltd

8/467/2A1, Nehru Street, 2nd Cross Road, 4th Lane, Arthi Nagar,  
Saraswathi Nagar, Thirumullaivaoyal, Chennai, Tamil Nadu 600062

+91 44 3559 4099 | +91 93811 27879

sales@itplmedical.com | www.itplmedical.com

Manufactured by

 **Huons  
Meditech**  
Korea



## FAST AND EFFECTIVE CYCLE TIME

- All Mode : 20 Minutes
- Quick Mode : 13 Minutes



## SELF DIAGNOSTIC SYSTEM

- Easily and clearly check the error code displayed on the device



## SELF-DISINFECTION MODE

- Completely disinfect tubes and chamber



## PERFECT CONTROL AND MANAGEMENT SYSTEM

- Fully customizable cycle options
- Cycle data recorded on print outs



## PROVIDING UNRIVALED RELIABILITY WITH TECHNOLOGICAL EFFICIENCY

- Creating a new standard on safe endoscope re-processing methods



## IMPROVED SAFETY

- The premium acrylic covers, specially designed to provide exceptional protection against both corrosion and breakage
- Soft closing chamber cover
- Cover sensor pauses cycle when the cover is open
- Equipped with powerful air-guns

### Product Information

Model	HUEN DR 02
Endoscope Quantity	2 Scope
Washing Method	3600 ROTATIONAL SPRAY
Chamber Capacity	16L
Disinfectant Capacity	18L
Detergent Tank Capacity	1L
Alcohol Tank Capacity	1L
Size (W*D*H)	696 x 656 x 1065 mm
Input Power	MAX 1,500W
Power Rating	AC220-230V (50/60 Hz)
Weight	81 kg
Air Gun & Printer	Standard

### Operating Procedure

#### Water + Detergent

1 min 30 Sec

#### Disinfection

5 min

#### Rinse

\*2 min

#### Alcohol

10 sec

#### Air

30 sec

#### Washing

00:00~90:50 (min;sec)

#### Disinfection

00:00~90:50 (min;sec)

#### Rinse

06:00~90:50 (min;sec)

\*Rinse time is fixed for safety

### Program(Automatic)

#### ALL MODE (20 min)

Wash - Detergent - Disinfect - Rinse - Alcohol - Dry

#### QUICK MODE (13 min)

Disinfect - Rinse - Alcohol - Dry