

Scope Buddy



Marketed & Promoted by

Innovative Therapeutics Pvt Ltd

8/467/2A1, Nehru Street, 2nd Cross Road, 4th Lane, Arthi Nagar,
Saraswathi Nagar, Thirumullaivaoyal, Chennai, Tamil Nadu 600062

+91 44 3559 4099 | +91 93811 27879

sales@itplmedical.com | www.itplmedical.com

Manufactured by

 **STERIS**

Product Overview

- Manual cleaning is essential in removing bioburden and preparing the endoscope for high-level disinfection. For optimal results,
- Cantel's manual cleaning products and detergents are designed to be a critical first step to efficiently and effectively eliminating bacteria on endoscopes prior to reprocessing.



Technical Specification

Unit Dimensions	9.5 H x 6 W x 8 D (inches) (24.1 H x 15.2 W x 20.3 D (cm))
Unit Weight	10 lbs (4.5 kg)
Electrical Requirements	110 - 230 V, 50/60 Hz
Electrical Safety Certifications	UL Standard 61010-1 CANADA/USA C22.2 No. 61010-1
Fluid Delivery Rate	1 liter/minute
Flush Time	0 - 99 minutes
Timer Control	Digital display with touchpad controls
System Pressure	Built-in safety bypass
Cycle Completion Indicators	Audible and visual alerts
Unit Placement	Counter-top or wall-mounting (bracket included)
Warranty	One year